
 **MONARK**
EXERCISE AB

MANUAL

Rehab Trainer
881E

English

Svenska

English

Congratulations on your decision to purchase a product manufactured by Monark in Sweden. Monark has been the world's leading manufacturer of high quality ergometers and exercise cycles for more than 40 years.

CONTENTS

Assembly instructions	4
Calibration	6
Operation instructions	7
How to replace batteries	11
Comparison of arm work to leg work	12
References	12

ASSEMBLY INSTRUCTIONS

The parts shown below are not assembled.


Fig 1

ASSEMBLY INSTRUCTIONS *cont.*

Assemble supporting tube with two screws M6 x 27 and two spring washers. See fig 2.

Pull out the pin and assemble the transport handle. Put the transport handle in its highest position when the Rehab Trainer is moved on the transport wheels. See fig 3.

Assemble crank arm marked R (Right) including two star locks and bolt M8 x 40 on the right hand side. Crank arm marked L (Left) including lock head and bolt M8 x 40 to be assembled on the left hand side.

Pedals

Assemble the pedal marked R (Right) on the right hand side. The pedal axle has a right hand thread and must be threaded onto the crank clockwise. The pedal marked L (Left) to be assembled on the left hand side. The pedal axle marked L has a left hand thread and must be threaded onto the crank counter clockwise.

Crank handle

The crank handle marked R (Right) is to be assembled on the right hand side with the hand guard placed on the axle. The axle has a right hand thread and must be threaded onto the crank clockwise. The crank handle marked L (Left) is to be assembled on the left hand side with the hand guard placed on the axle. The axle has a left hand thread and must be threaded onto the crank counter clockwise.

NOTE! The pedals or the crank handles can be assembled in three different positions on the crank arms. See page 6, fig 7.


Fig 2


Fig 3


Fig 4

CALIBRATION

Calibration is done at the factory. If you for some reason want to check the calibration, do as follows:

Tighten the brake belt through turning the control knob so that the pointer goes up to about 25-50 Watt at rotation of the crank.

Put the Rehab Trainer at the edge of a table as shown in fig 5.

Loosen the screws for the clamping plate, so that the end of the brake belt is loose. See fig 5.

Fasten a 2 kg weight (our Article No 9000-212) to the brake belt. See fig 5. This weight should now be read on the scale at 2 kp. See fig 6.

Should there be a deviation between the position of the pointer and the 2 kp mark on the scale, adjust the loading spring by turning the adjusting screw until the correct position is obtained. See fig 5 and 6. Turn clockwise if the pointer is placed too high and turn counter clockwise if the pointer is placed too low.

Remove the weight and fasten the brake belt with the screw for the clamping plate so that the pointer does not go below the zero ("0") mark on the scale.

The height level of the cranks is adjustable when loosening the nuts according to fig 7. Set the crank arms at the desired position and fasten with the nuts.

If the cranks are not in line, adjust these. Set first the left hand crank straight upwards. Then loosen the right hand crank by undoing the screw as per fig 7. Adjust the position of the right hand crank so that it is in line with left hand crank and fasten the screw.


Fig 5


Fig 6


Fig 7

OPERATION INSTRUCTIONS

Monark Rehab Trainer model 881E is an arm and leg ergometer provided with a belt brake. The power can be read in Watts at 50 pedal revolutions per minute. Moreover, the work can be read in kilopondmeters (kpm). The Rehab Trainer is also equipped with an electronic meter, showing pedal revolutions per minute, the total pedal revolutions and time function.

When cycling the test person puts supplies the flywheel into a certain kinetic energy. This is braked by means of a brake belt which runs around the bigger part of the brake surface of the flywheel.

The workload is changed either by using another pedalling speed or by increasing or decreasing the tension of the brake belt against the flywheel by means of the work load control knob.


Fig 12


Fig 8 Arm exercise in a sitting position with the ergometer placed on a table


Fig 9 Arm exercise with the ergometer hanging on wall bars.


Fig 10 Leg exercise in a sitting position with the ergometer placed on the floor.


Fig 11 Leg exercise in a lying position at the same level as the ergometer.

OPERATION INSTRUCTIONS *cont.*

The total work can be read in kilopondmeter. In our example we have chosen 50 watt as work load. When 250 revolutions are finished, we can read the result that the work is about 1500 kilopondmeter.

PEDAL REVOLUTIONS												
<i>Chosen work load</i>												<i>Total pedal revolutions</i>
	25	50	75	100	125	150	175	200	225	250	275	300
	kpm											
25	75	150	225	300	375	450	525	600	675	750	825	900
50	150	300	450	600	750	900	1050	1200	1350	1500	1650	1800
75	225	445	675	900	1125	1350	1575	1800	2025	2250	2475	2700
100	300	600	900	1200	1500	1800	2100	2400	2700	3000	3300	3600

Fig 13

Meter functions

The buzzer gives an audible signal after half the exercise time and when there are 30 seconds left. It gives repeated signals for 20 seconds after the end of the pre-set time.

Start by pressing the button for function. Press as many times as is necessary to get the marker on "REV" according to fig 15. To reset, press the reset button so that display according to fig 15 is shown.


Fig 14

OPERATION INSTRUCTIONS *cont.*

Press the function button once so that the marker stands at "TIMER" according to fig 16. If the display is not reset, press the reset button so that a display like fig 17 is shown. After 5 seconds the display changes automatically to fig 16.

Let's imagine that you will exercise for 1 hour and 15 minutes. Press the buttons for hours and minutes, so that a display like fig 18 is shown. The countdown starts 5 seconds after the time has been entered.

When setting time exceeding one hour, the meter shows time in hours and minutes but not in seconds, which are then existing only in the computer. This means that when the countdown starts after about 5 seconds, the display shows 1.14 at the same time as the computer invisibly, counts the remaining seconds up to 1 hour and 14 minutes.

At time below one hour both minutes and seconds are shown on the display.

Press the function button once, so that the display shows "SCAN". See fig 19.

The displaying will now automatically change between the function at intervals of 3.5 seconds.

Scan cycles through the functions as shown in fig. 20.


Fig 15


Fig 16


Fig 17


Fig 18


Fig 19


Fig 20

OPERATION INSTRUCTIONS *cont.*

When the marker stand at "RPM" pedal revolutions per minute are shown. Fig 21 shows 50 pedal revolutions per minute.

When the marker stands at "REV" the total number of revolutions achieved during the exercise is shown. Fig 22 shows a total of 500 pedal revolutions.

5 minutes after the cycling is finished, the meter is automatically shut off.

After the exercise is finished, the registered total number of pedal revolutions are kept in the meter.

If the exercise is stopped more than 5 minutes before the exercise time setting, the remaining time – when the meter is automatically shut off – will still be there, when the meter is started over again.

As to resetting, please see figures 15, 16 and 17 at page 9.

WARNING! Do not expose the electronic meter to direct sunlight or extremely high temperatures. Do not use any dissolvents when cleaning but a dry cloth.

Handling

Work is started with a slack brake belt. Thereafter, the belt should be stretched with the tensioning knob until the required workload is obtained. A metronome will assist the subject in maintaining the proper rotational speed. As the belt and wheel get warmed up, the friction will change, necessitating readjustment, especially if the trainer has been unused for any length of time.

Published test protocols for the arm ergometer have been limited. One of the few in print by Jones and Campbell, recommended increments of 50 kpm/min (approximately 8 Watts) for their Stage 1 study. (5)


Fig 21


Fig 22

Table 1. Oxygen Consumption (L/min) at Workloads Below the Anaerobic Threshold:

Work (kpm/min)	(Approx Watts)	V02	Work (kpm/min)	(Approx Watts)	V02
150	24	.6	400	64	1.1
200	32	.7	450	72	1.2
250	40	.8	500	80	1.3
300	48	.9	550	86	1.4
350	56	1.0	600	96	1.5

HOW TO REPLACE BATTERIES

Remove the sliding lid and replace the batteries. Put the sliding cover back in its position.


Fig 23


Fig 24


Fig 25

COMPARISON OF ARM WORK TO LEG WORK

When comparing arm ergometry or arm exercise to leg ergometry, there are several physiologic considerations. The maximum work that a muscle is able to perform is related to its size. Obviously, the muscle mass of the arms is considerably smaller than the legs and therefore the maximum work and maximum oxygen consumption of the legs will be greater. Hershfield and others added that the mechanical differences in the shoulder and hip joints increase metabolic demands when arm work is performed. (1)

Reybrouck and others found that at workloads below the "anaerobic threshold", when adequate oxygen is supplied to meet all the muscle's requirements, the oxygen consumption from arm work matches closely the oxygen consumption during equivalent leg work. (2) They noted that heart rate and ventilation were higher during arm work at the same oxygen consumption in leg ergometry.

In rehabilitation or evaluation of disability, the expected sustained work capability of an individual is approximately 40% of maximum work and below the anaerobic threshold. (3) To use the equivalent work tables which have been published (reference 4), divide the V_{O_2} at the subject's workload by the patient's weight in kilograms. Divide this V_{O_2} ml/kg by 3.5. This will give the equivalent METs and then evaluate the tasks of which they are capable of performing at that MET level.

REFERENCES

1. Hershfield S, et al. Relative effects on the heart by muscular work in the upper and lower extremities. *Arch Phys Med Rehab*. Pages 249-257, May 1968.
2. Reybrouck T, et al. Limitations to maximum oxygen uptake in arm, leg and combined arm-leg ergometry. *J Appl Physiol* 38:774-669, 1975.
3. Becklage MR. Organic or functional impairment. *Am Rev Resp Dis* 129:S96-100, 1984.
4. Physician's Handbook for Evaluation and Physical Fitness. Tennessee Heart Association, Physical Exercise Committee, 1972, pgs 44-45.
5. Jones NL, Campbell EJM. Clinical exercise testing. WB Saunders Publishing, Phila, PA. Page 86, 1975.

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