

Since 1946, Whitehall's reputation for quality continues to be trusted in health care facilities around the world. All products are thoughtfully engineered to provide the most comfortable and effective treatment for the patient and reliable convenience to the caregiver.



Professional Health Care and Therapy Equipment

WHIRLPOOLS

Whitehall Whirlpools are flexible in their application, durable in their construction, and effective in their design. Temperatures as well as pressures exerted from water carry impulses through the nerves into the body, which is beneficial to blood flow, circulation, and pain sensitivity.

DAKON



**Whitehall
WHIRLPOOLS**



Stationary units are faster to fill and easily drained for high-traffic therapy centers where time is a factor

Mobile whirlpools allow for convenient filling, use and storage




Dual-scale (Fahrenheit and Celsius) thermometer with adjustable calibration

Fully-covered bottom prevents against bacteria buildup



WHIRLPOOLS

Our whirlpools are:

- Durable, 16-gauge type 304 **stainless steel construction**
- 80% **recycled** material 
- Backed by a **one-year warranty**

Powerful $\frac{1}{2}$ horsepower turbine **motor**, MET compliant, is standard on all whirlpool models

Rolled stainless steel tank **rims** for added comfort



Anti-rattle retainer for thermometer

Seamless tanks are ground and polished to a **satın finish**

Heavy-duty **galvanized steel** and **rubber casters**

